

|                                                                |                          |              |
|----------------------------------------------------------------|--------------------------|--------------|
| <b>Caesiumchlorid p.A.</b><br>Cesium chloride analytical grade | <b>Kat.Nr./Cat.No. :</b> | <b>15554</b> |
|----------------------------------------------------------------|--------------------------|--------------|

| <b>Parameter</b><br>parameter               | <b>Methode</b><br>method | <b>Spezifikation</b><br>specification       |
|---------------------------------------------|--------------------------|---------------------------------------------|
| <b>Summenformel</b><br>formula              |                          | CsCl                                        |
| <b>Molekulargewicht</b><br>molecular weight |                          | 168.36                                      |
| <b>Aussehen</b><br>appearance               |                          | weißes krist. Pulver<br>white cryst. powder |
| <b>Gehalt (%)</b><br>assay                  |                          | min. 99.9                                   |
| <b>A 1cm/ 260nm</b>                         | 50% in H <sub>2</sub> O  | max. 0.02                                   |
| <b>A 1cm/ 280nm</b>                         | 50% in H <sub>2</sub> O  | max. 0.02                                   |
| <b>Schwermetalle (ppm)</b><br>heavy metals  | als/ as Pb               | max. 1                                      |
| <b>Lagerung (°C)</b><br>storage             |                          | +15 bis/ to +30                             |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
We do not guarantee that the product can be used for a special application.