

Carbenicillin-Na₂-salz reinst Carbenicillin-Na ₂ -salt research grade	Kat.Nr./Cat.No. : 15875
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₇ H ₁₆ N ₂ O ₆ S·Na ₂
Molekulargewicht molecular weight		422.37
Aussehen appearance		weißes bis fast weißes Pulver white to off-white powder
Identität identity	IR	entspricht corresponds
Aktivität (µg/ mg) activity	bez. auf getr. Subst. calc. on dried subst.	min. 770
pH	1 % in H ₂ O	6.0 – 7.5
Wasser (%) water	KF	max. 5.0
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.