

<b>Carmin</b> reinst Carmine research grade	<b>Kat.Nr./Cat.No. :</b>	<b>16180</b>
--	--------------------------	--------------

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Aussehen</b> appearance		grau-violettes Pulver grey-violet powder
<b><math>\lambda</math> max (nm)</b>	0.001 % in 0.1N NaOH	528 – 536
<b>A 1cm/ <math>\lambda</math>max</b>	0.001 % in 0.1N NaOH	min. 0.045
<b>Wasser (%)</b> water	KF	max. 16.0
<b>Löslichkeit</b> solubility		schwer löslich in Wasser slightly soluble in water löslich in Ammoniak soluble in ammonia
<b>TLC</b>	Kieselgel/silicagel 60 F <sub>254</sub> Aceton/ acetone 10 Ammoniak/ ammonia (25 %) 10	einheitlich one spot
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
We do not guarantee that the product can be used for a special application.