

Cefotaxim·Na-Salz reinst Cefotaxime·sodium-salt research grade	Kat.Nr./Cat.No. : 16382
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{16}H_{16}N_5O_7S_2 \cdot Na$
Molekulargewicht molecular weight		477.5
Aussehen appearance		weißes Pulver white powder
Gehalt (%) content		1g Cefotaxim·Na entspricht 954.2 mg Cefotaxim 1g Cefotaxime·sodium-salt contains 954.2 mg Cefotaxime
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.