

Cholsäure-Na-salz p.A. Cholic acid-Na-salt analytical grade	Kat.Nr./Cat.No. :	17126
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{24}H_{39}O_5 \cdot Na$
Molekulargewicht molecular weight		430.5
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration	98.0 – 101.0
$[\alpha]_D/20^\circ C$	10% in H ₂ O	+26° bis/ to +28°
Trocknungsverlust (%) loss on drying		max. 5.0
pH	10% in H ₂ O	8.0 – 9.5
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 20
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.