

Cholin-chlorid reinst Choline-chloride research grade	Kat.Nr./Cat.No. :	17130
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_5H_{14}ON^+ \cdot Cl^-$
Molekulargewicht molecular weight		139.6
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	Titration	98.0 –100.5
Identität identity	IR	entspricht corresponds
Trocknungsverlust (%) loss on drying	120 °C	max. 1.5
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 10
Glührückstand (%) residue on ignition		max. 0.1
Ursprung origin		nicht tierisch non animal
Lagerung (°C) storage		+15 bis/ to +30

Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen des DAC.
The physical/ chemical data corresponds to the requirements of DAC.

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--