
Collagenase N Reference Identity from <i>Clostridium histolyticum</i> ; EC 3.4.24.3 lyophil.	Cat. No. :	17464
--	-------------------	--------------

Parameter	Specification
------------------	----------------------

Characters

Constitution	Fine powder
Colour	Light brown

Identification with SDS-PAGE

Pattern of the main bands	Must comply
Collagenase AZ	Must comply
Clostridiopeptidase A	Must comply
Clostripain heavy chain	Must comply
Neutral protease	Must comply
Clostripain light chain	Must comply

Tests

Ash	• 3.0 %
Sulphur	• 2.0 %

Assay

Clostridiopeptidase A activity	Must comply
--------------------------------	-------------

Use: not for use in human beings
Store at -15 to -25 °C

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.