

Coomassie® Brilliant Blau R 250 rein Coomassie® Brilliant Blue R 250 pure	Kat.Nr./Cat.No. :	17525
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄₅ H ₄₄ N ₃ O ₇ S ₂ Na
Molekulargewicht molecular weight		826.0
Aussehen appearance		schwarz-rotes Pulver black to red-colored powder
Gehalt (%) assay	UV	min. 75.0
Wasser (%) water	KF	max. 10.0
TLC		entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.