

D-Cylcoserin reinst D-Cylcoserine research grade	Kat.Nr./Cat.No. :	17685
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₃ H ₆ N ₂ O ₂
Molekulargewicht molecular weight		102.1
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay		min. 98.0
λ_{max}	0.001% in H ₂ O	224 - 228
A 1cm/λ_{max}	1% in H ₂ O	min. 300
[α]_D/20°C	1% in H ₂ O	+110° bis/ to +118°
pH	10% in H ₂ O	5.7 – 6.3
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.