

---

<b>Cytochrom C aus Pferdeherz lyophil.</b> Cytochrome C from equine heart lyophil.	<b>Kat.Nr./Cat.No. :</b>	<b>18022</b>
---	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Molekulargewicht</b> molecular weight		ca. 12 400
<b>Aussehen</b> appearance		oxidrotes Lyophil. oxid-red lyophil.
<b>Löslichkeit</b> solubility	3 mg/ ml in H <sub>2</sub> O	klare dunkelrote Lösung clear dark-red solution
<b>Trocknungsverlust (%)</b> loss on drying		max. 10.0
<b>Gehalt reduziert (%)</b> assay reduced	photometrisch photometric	90 – 100
<b>Lagerung (°C)</b> storage		-15 bis/ to -25

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.