

Desoxycholsäure·Na-salz rein Deoxycholic acid·Na-salt pure	Kat.Nr./Cat.No. :	18330
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{24}H_{39}O_4 \cdot Na$
Molekulargewicht molecular weight		414.6
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay		min. 99.0
A 1cm/ 230nm	1% in H ₂ O	max. 0.5
A 1cm/ 260nm	1% in H ₂ O	max. 0.2
A 1cm/ 280nm	1% in H ₂ O	max. 0.1
pH	1% in H ₂ O	7.5 – 9.5
Trocknungsverlust (%) loss on drying		max. 5.0
Natriumcholat (%) sodium cholate		max. 2.0
Schwermetalle (ppm) heavy metals		max. 50
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.