

Cyanase(TM) Nuclease Lösung in 50 % Glycerin
Cyanase(TM) Nuclease solution in 50 % glycerol

Cat. No. :

18542

Parameter parameter	Methode method	Spezifikation specification
Reinheit purity	SDS – PAGE	entspricht corresponds
Enzymkonzentration (U/ µl) concentration		50
Protease Aktivität protease activity	Fluoreszenz-basierter Test fluorescent based assay	nicht nachweisbar not detectable
Lagerung (°C) storage		-15 bis/ to -25

Unit definition:

One unit is the amount of enzyme that degrades 3 µg of full length λ DNA completely at 37 °C in one minute.
Reaction conditions are 50 mM Tris pH 8.0, 6 mM MnSO₄.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.