

4',6-Diamidino-2-phenylindol-2HCl p.A.	Kat.Nr./Cat.No. :	18860
4',6-Diamidino-2-phenylindole-2HCl analytical grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{16}H_{15}N_5 \cdot 2HCl$
Molekulargewicht molecular weight		350.25
Aussehen appearance		gelbes Pulver yellow powder
Reinheit (%) purity		min. 98
Löslichkeit solubility	2.5 % in H ₂ O	klare, gelbe Lösung clear, yellow solution
λ_{max} (nm)	0.001 % in H ₂ O	340 – 348
A 1 cm/ λ_{max}	0.001 % in H ₂ O	min. 0.7
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.