

2,6-Dichlorphenolindophenol-Na p.A.
Tillmans Reagens analytical grade

Kat.Nr./Cat.No. : 19354

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{12}H_6Cl_2NO_2 \cdot Na$
Molekulargewicht molecular weight		290.09
Aussehen appearance		grau-schwarzes Pulver mit Grünstich greyish black powder with a greenish cast
λ_{max}	0.001% in H ₂ O	600 – 608nm
A 1cm/ λ_{max}	0.001% in H ₂ O	min. 0.4
Trocknungsverlust (%) loss on drying		max. 12.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.