

<b>Digitonin p.A.</b> (Digitin) analytical grade	<b>Kat.Nr./Cat.No. :</b>	<b>19550</b>
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>56</sub> H <sub>92</sub> O <sub>29</sub>
<b>Molekulargewicht</b> molecular weight		1 229,34
<b>Aussehen</b> appearance		weißes Pulver white powder
<b>Löslichkeit</b> solubility	USP	entspricht corresponds
<b>Trocknungsverlust (%)</b> loss on drying		max. 6.0
<b>[α]<sub>D</sub>/20°C</b>	10% in 75% HAc	-47.0° bis/ to -49.0°
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen der USP.**  
The physical/ chemical data corresponds to the requirements of USP.

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--