

Dithioerythrit reinst DTE research grade	Kat.Nr./Cat.No. :	20697
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₄ H ₁₀ O ₂ S ₂
Molekulargewicht molecular weight		154.24
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration	min. 99.0
Schmelzpunkt (°C) melting point		82 – 86
A 1cm/ 283nm	0.02 M in H ₂ O	max. 0.1
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.