

Dodecylsulfat·Na-Salz 2x krist. p. A.	Kat.Nr./Cat.No. :	20763
Dodecylsulfate·Na-salt 2x cryst. analytical grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{12}H_{25}O_4S \cdot Na$
Molekulargewicht molecular weight		288.38
Aussehen appearance		weißes krist. Schuppen oder Pulver white cryst. flakes or powder
Gehalt (%) assay	GC	
	C ₁₀ -Sulfat / sulfate	< 1
	C ₁₂ -Sulfat / sulfate	min. 99.0
	C ₁₄ -Sulfat / sulfate	< 1
Löslichkeit solubility	10% in H ₂ O	klar, max. 20 Hazen clear, max. 20 Hazen
A 1cm/ 260 nm	10% in H ₂ O	max. 0.2
A 1cm/ 280 nm	10% in H ₂ O	max. 0.08
Sulfatasche (%) sulfated ash		21.5 – 27.5
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.