

**Dodecylsulfat·Na salz, Pellets reinst** **Kat.Nr./Cat.No. : 20765**  
Dodecylsulfate·Na salt , in pellets research grade

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>12</sub> H <sub>25</sub> O <sub>4</sub> SNa
<b>Molekulargewicht</b> molecular weight		288.38
<b>Aussehen</b> appearance		weiße bis fast weiße Pellets white to off-white pellets
<b>Gehalt (%)</b> assay	GC	> 99.0
	C <sub>10</sub> -Sulfat / sulfate	< 1
	C <sub>12</sub> -Sulfat / sulfate	min. 99.0
	C <sub>14</sub> -Sulfat / sulfate	< 1
<b>Löslichkeit</b> solubility	10% in H <sub>2</sub> O	klar, max. 50 Hazen clear, max. 50 Hazen
<b>A 1cm/ 260 nm</b>	10% in H <sub>2</sub> O	max. 1.5
<b>A 1cm/ 280 nm</b>	10% in H <sub>2</sub> O	max. 1.0
<b>Sulfatasche (%)</b> sulfated ash		21.5 – 27.5
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.