

<b>Dodecylsulfat·Na-Salz 2x krist.</b> Dodecylsulfate·Na-salt 2x cryst. premium grade	<b>Kat.Nr./Cat.No. :</b>	<b>20770</b>
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>12</sub> H <sub>25</sub> O <sub>4</sub> S·Na
<b>Molekulargewicht</b> molecular weight		288.38
<b>Aussehen</b> appearance		weißes Pulver white powder
<b>Gehalt (%)</b> assay	GC	min. 99.0
<b>Löslichkeit</b> solubility	10% in H <sub>2</sub> O	klar, max. 5 Hazen clear, max. 5 Hazen
<b>A 1cm/ 260 nm</b>	10% in H <sub>2</sub> O	max. 0.15
<b>A 1cm/ 280 nm</b>	10% in H <sub>2</sub> O	max. 0.05
<b>Sulfatasche (%)</b> sulfated ash		21.5 - 27.5
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--