

Eastman Instant Adhesive 910

Kat.Nr./Cat.No. :

20832

| Parameter | | Specification |
|--|-----------|------------------------|
| Specific gravity (g cm ⁻³) | | 1.02 – 1.08 |
| Flash point (°C) | | 82 |
| Flexural modulus (MNm ⁻²) | D970 | 3400 |
| Water absorption (%) | 25°C | 0.6 |
| Vapour pressure (mm Hg) | liquid | 2.2 |
| Dissipation factor (HZ) | D150 | 0.02 x 10 ⁶ |
| Dielectric constant (HZ) | D150 | 3.34 x 10 ⁶ |
| Volume resistivity (ohm cm) | D257 | 10 ¹⁴ |
| Thermal conductivity (Cal/cm sec °C) | | 7.2 x 10 ⁻⁴ |
| Boiling point (°C) | 2.0 mm Hg | 48 – 54 |
| Storage (°C) | | +2 to +8 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.