

Eosin Y·Na-Salz reinst	Kat.Nr./Cat.No. :	21005
Eosin Y·Na-salt research grade, C.I. 45380		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{20}H_6Br_4O_5 \cdot Na_2$
Molekulargewicht molecular weight		691.9
Aussehen appearance		rotbraunes Pulver reddish brown powder
λ max (nm)	0.001 % in H ₂ O	516 ± 4
A 1cm/ λ max	0.001 % in H ₂ O	min.1.0
ϵ λ max		min. 69 000
Wasser (%) water	KF	max. 15.0
TLC	Kieselgel/silicagel 60 F ₂₅₄ CH ₂ Cl ₂ 7 Methanol 3 H ₂ O 0.5	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.