
Embedding Medium ERL-4221D pract. Cycloaliphatic Epoxy Resin 3,4-Epoxycyclohexylmethyl-3,4-epoxycyclohexylcarboxylat	Kat.Nr./Cat.No. : 21041
---	--------------------------------

Parameter parameter	Spezifikation specification
Summenformel formula	$C_{14}H_{20}O_4$
Molekulargewicht molecular weight	252.3
Aussehen appearance	klare, farblose Flüssigkeit clear, colorless liquid
Farbe (Pt-Co) color	max. 50
Dichte (g/ml) density	1.159 – 1.174
Viskosität (mPa*s at 25°C) viscosity	220 – 250
Epoxy Äquivalentgewicht epoxy equivalent weight	126 – 135
Lagerung (°C) storage	+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.