

<b>Ethidiumbromid</b> reinst Ethidium bromide research grade	<b>Kat.Nr./Cat.No. :</b>	<b>21238</b>
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		$C_{21}H_{20}N_3^+ \cdot Br^-$
<b>Molekulargewicht</b> molecular weight		394.33
<b>Aussehen</b> appearance		rotes Pulver red powder
<b>Gehalt (%)</b> assay	Titration	min. 98.0
<b>Trocknungsverlust (%)</b> loss on drying		max. 5.0
$\lambda_{max}$ (nm)	0.001 % in Methanol	523 – 528
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.