

Fast Green FCF research grade	Kat.Nr./Cat.No. : 21295 Lot/Contr.No. :
---	--

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{37}H_{34}N_2O_{10}S_3 \cdot Na_2$
Molekulargewicht molecular weight		808.9
Aussehen appearance		weinrotes Pulver purple powder
λ_{max} (nm)	0.0005% in H ₂ O	621 - 627
A 1cm / λ_{max}	0.0005% in H ₂ O	min. 0.6
Wasser (%) water	KF	max. 20.0
TLC	Kieselgel/silicagel 60 F ₂₅₄ n-Butanol/n-butanol 5 Ethanol/ethanol 2 Eisessig/glacial acetic acid 1 Wasser/water 2	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.