

Fast Red TR-Salz Fast Red TR-salt	Kat.Nr./Cat.No.:	21317
---	-------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_7H_6ClN_2 \cdot C_{10}H_7O_6S_2$
Molekulargewicht molecular weight		440.9
Aussehen appearance		weißes bis leicht gelbliches Pulver white to light yellow powder
λ_{max} (nm)	0.001 % in H ₂ O	286 - 290
A 1cm/ λ_{max}	0.001 % in H ₂ O	min. 0.25
Wasser (%) water	KF	max. 10.0
TLC	Kieselgel/silicagel 60 F ₂₅₄ n-Butanol Eisessig/acetic acid Wasser/water	entspricht corresponds 40 20 20
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.