

Ferrozine[®] p.A. Ferrozine [®] analytical grade	Kat.Nr./Cat.No. :	21326
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{20}H_{13}N_4O_6S_2 \cdot Na$
Molekulargewicht molecular weight		492.5
Aussehen appearance		hellgelbes Pulver light yellow powder
Gehalt (%) assay	Titration	min. 98.0
λ_{max} (nm)	0.001% in H ₂ O	284 ± 2
$\epsilon\lambda_{max}$	in H ₂ O	min. 20 000
Wasser (%) water	KF	max. 10.0
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.