
Fibronectin (human) lyophil., analytical grade, sterile	Kat.Nr./Cat.No. :	21370
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Molekulargewicht molecular weight		ca. 450 000
Aussehen appearance		weißes Lyophilisat white lyophilisate
Endotoxine (ng/mg protein)	LAL	
Reinheit (%) purity	SDS-Elektrophorese SDS-electrophoresis	> 90
Sterilität sterility		negativ negative
HBsAg , anti-HIV 1 und 2 HBsAg , anti-HIV 1 and 2		negativ negative
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.