

Flavin-adeninucleotid-Na₂-salz rein FAD·Na ₂ pure	Kat.Nr./Cat.No. :	21490
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₇ H ₃₁ N ₉ O ₁₅ P ₂ ·Na ₂ ·2H ₂ O
Molekulargewicht molecular weight		865.56
Aussehen appearance		helloranges Pulver light orange powder
Gehalt (%) assay	aus/ from ε 450nm (bez auf getr. Subst./ rel. on dried subst.)	min. 95.0
pH	1% in H ₂ O	6.2 ± 0.5
Identität identity	IR	entspricht corresponds
Wasser (%) water	KF	max. 10.0
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.