

Fluorescein diacetat 5xkrist. reinst Fluorescein diacetat 5xcryst. research grade	Kat.Nr./Cat.No. :	21575
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₄ H ₁₆ O ₇
Molekulargewicht molecular weight		416.4
Aussehen appearance		leicht gelbes Pulver pale yellow powder
Schmelzpunkt (°C) melting point		200 – 205
Identität identity	IR	entspricht corresponds
Dünnschichtchromatographie thin-layer chromatography	TLC	entspricht corresponds
Sulfatasche (%) sulfated ash		max. 0.1
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.