

D-Fructose-1,6-biphosphat-Na₃-salz , reinst D-Fructose-1,6-biphosphate-Na ₃ -salt (Harden-Young-Ester), research grade	Kat.Nr./Cat.No. :	21860
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_6H_{11}O_{12}P_2Na_3 \cdot 8 H_2O$
Molekulargewicht molecular weight		550.19
Aussehen appearance		reinweißes Pulver off-white powder
Gehalt (%) assay	Titration	min. 98.0
Identität identity	IR	entspricht corresponds
Wasser (%) water	KF	ca. 25.0
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.