

|                                                   |                          |              |
|---------------------------------------------------|--------------------------|--------------|
| <b>Fuchsin basisch</b> rein<br>Fuchsin basic pure | <b>Kat.Nr./Cat.No. :</b> | <b>21916</b> |
|---------------------------------------------------|--------------------------|--------------|

| Parameter<br>parameter                   | Methode<br>method                                                                                          | Spezifikation<br>specification              |
|------------------------------------------|------------------------------------------------------------------------------------------------------------|---------------------------------------------|
| <b>Aussehen</b><br>appearance            |                                                                                                            | grünes krist. Pulver<br>green cryst. powder |
| <b><math>\lambda_{\max}</math> (nm)</b>  | 0.0001% in 50% EtOH                                                                                        | 552 $\pm$ 4                                 |
| <b>A 1cm/<math>\lambda_{\max}</math></b> | 0.0001% in 50% EtOH                                                                                        | min. 0.15                                   |
| <b>Wasser (%)</b><br>water               | KF                                                                                                         | max. 15.0                                   |
| <b>TLC</b>                               | Nano-Sil C <sub>18</sub> -100 UV <sub>254</sub><br>MeOH 3<br>NH <sub>3</sub> (25%) 1<br>H <sub>2</sub> O 2 | entspricht<br>corresponds                   |
| <b>Lagerung (°C)</b><br>storage          |                                                                                                            | +15 bis/ to +30                             |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
We do not guarantee that the product can be used for a special application.