

D-Galactosamin·HCl reinst	Kat.Nr./Cat.No. :	22009
D-Galactosamine·HCl research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_6H_{13}NO_5 \cdot HCl$
Molekulargewicht molecular weight		215.6
Aussehen appearance		weißes krist. Pulver white cryst. powder
Gehalt (%) assay	Titration	min. 98.0
Identität identity	IR	entspricht corresponds
Trocknungsverlust (%) loss on drying		max. 1.0
Schmelzpunkt (°C) melting point		176 – 181
[α]D/20°C	1% in H ₂ O nach / after 24h	+93.0° bis/ to +99.0°
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.