

γ -Globulin bovine
lyophil.

Cat.No. : 22550

Parameter	Method	Specification
Protein (%) (protein factor 6.25)	Nitrogen Analyzer dry weight basis	96 – 100
Purity (%) (γ -globulins, 7% solution)	Cellulose Acetate Electrophoresis	96 – 100
pH	7% solution	6.8 – 7.2
Sodium (mg/g)	Flame Photometer	0 – 10.0
Chloride (mg/g)	Coulometric Titrator	0 – 24.0
Moisture (%)	KF	0 – 5.0
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.