

---

<b><math>\alpha</math>-D-Glucose</b> p.A., wasserfrei Dextrose, analytical grade, anhydrous	<b>Kat.Nr./Cat.No. :</b>	<b>22700</b>
--	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>6</sub> H <sub>12</sub> O <sub>6</sub>
<b>Molekulargewicht</b> molecular weight		180.16
<b>Aussehen</b> appearance		weißes Pulver white powder
<b>Gehalt (%)</b> assay	Titration	min. 99.5
<b>Trocknungsverlust (%)</b> loss on drying		max. 0.5
<b>pH</b>	10% in H <sub>2</sub> O	5.0 – 7.0
<b>Sulfatasche (%)</b> sulfated ash		max. 0.05
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.