

D-Glucuronsäure·Na-salz reinst	Kat.Nr./Cat.No. :	22890
D-Glucuronic acid·Na-salt research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_6H_9O_7 \cdot Na \cdot H_2O$
Molekulargewicht molecular weight		234.1
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	Titration	min. 99.0
Löslichkeit solubility	5% in H ₂ O	klar und farblos clear and colorless
Wasser (%) water	KF	max. 10.0
$[\alpha]_D^{20}$	2% in H ₂ O nach/ after 2h	+21.0° ± 1°
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis / to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--