

---

<b>Glutaraldehyd 25% in Wasser</b> Standardqualität Glutaraldehyde 25% in water standard grade	<b>Kat.Nr./Cat.No. :</b>	<b>23115</b>
---	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>5</sub> H <sub>8</sub> O <sub>2</sub>
<b>Molekulargewicht</b> molecular weight		100.13
<b>Aussehen</b> appearance		klare, farblose Lösung clear, colorless solution
<b>Gehalt (%)</b> assay	RI	25.0 – 26.0
<b>Brechungsindex</b> refractive index	20°C	1.3720 – 1.3735
<b>Purification index</b>	A 235 / A 280; 1% in H <sub>2</sub> O	max. 5.0
<b>pH</b>	25% in H <sub>2</sub> O	3.0 – 4.0
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.