

Glutaraldehyd 50% in Wasser für Elektronenmikroskopie Glutaraldehyd 50% solution in water for electron microscopy	Kat.Nr./Cat.No. : 23116
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₅ H ₈ O ₂
Molekulargewicht molecular weight		100.13
Aussehen appearance		klare, farblose Lösung clear, colorless solution
Gehalt (%) assay	Titration	45.0 – 53.0
Brechungsindex refractive index	20°C	1.410 – 1.421
Purification index	A 235/ A280 1% in H ₂ O	
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.