

L-Glutathion oxidiert krist. reinst L-Glutathione (oxidized form) cryst. research grade	Kat.Nr./Cat.No. :	23130
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{20}H_{32}N_6O_{12}S_2$
Molekulargewicht molecular weight		612.64
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	HPLC	min. 98.0
Trocknungsverlust (%) loss on drying	3 h, 105 °C	max. 8.0
[α]_D/ 20 °C	4 % in H ₂ O	-106.0° bis/ to -96.0°
IR		entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.