

L-Glutathion reduziert krist. reinst L-Glutathione (reduced form) cryst. research grade	Kat.Nr./Cat.No. : 23150
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₇ N ₃ O ₆ S
Molekulargewicht molecular weight		307.3
Aussehen appearance		weißes Pulver white powder
Reinheit (%) purity	CE	min. 97.5
verwandte Substanzen (%) related substances	GSSG	max. 1.5
Schmelzpunkt (°C) melting point		182 – 192
Trocknungsverlust (%) loss on drying	3 h, 105 °C	max. 0.5
[α]_D/20°C	4 % in H ₂ O	-15.5° bis/ to -17.5°
Sulfatasche (%) sulfated ash		max. 0.1
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 10
Ursprung origin		nicht tierisch non animal
Lagerung (°C) / storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.