

Glycerin-2-phosphat Dinatriumsalz Pentahydrat rein 2-Glycerolphosphate-Na ₂ -salt Pentahydrate pure	Kat.Nr./Cat.No. :	23330
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₃ H ₇ O ₆ P·Na ₂ ·5H ₂ O
Molekulargewicht molecular weight		306.1
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	α + β-Glycero-PO ₄	min. 99.0
Wasser (%) water	KF	26 – 33
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.