
Guanosin-5'-triphosphat Na₂-salz p.A. **Kat.Nr./Cat.No. : 24321**
Guanosine-5'-triphosphate Na₂-salt analytical grade

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₄ N ₅ O ₁₄ P ₃ ·Na ₂ ·x 3H ₂ O
Molekulargewicht molecular weight		567.2 (wasserfrei/ anhydrous)
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	ε 252nm, pH 7.0 (bez. a. getr. Subst./ calc. on dried subst.)	min. 96.0
Wasser (%) water	KF	max. 10.0
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.