

Hämin reinst Hemin research grade	Kat.Nr./Cat.No. :	24410
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{34}H_{32}FeN_4O_4 \cdot Cl$
Molekulargewicht molecular weight		652.0
Aussehen appearance		schwarzblaues glänzendes Pulver dark-blue gleaming powder
Gehalt (%) assay	UV/575nm	min. 98.0
A 1cm/ 575 nm	0.004% in 0.1M NaOH cont. 2.5% Triton (Hämoglobin Reagenz hemoglobine reagent)	0.427 ± 0.004
Wasser (%) water	Trocknungsverlust loss on drying	max. 2.0
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.