

Hyaluronidase from ovine testes
lyophil. EC 3.2.1.35

Cat.No. : 25118

Parameter	Method	Specification
Molecular weight		ca. 55 000
Appearance		white lyophilisate
Solubility	1 mg/ml 0.02M Phosphate buffer pH 6.9 cont. 0.45% NaCl and 0.01% BSA	clear and colorless
Activity (U/mg)	Hyaluronidase	min. 1 000
Storage (°C)		-15 to -25

Unit definition

1 Unit is that amount of enzymatic activity which produces the same reduction in turbidity in a mixture of hyaluronic acid and albumin as one I.U. (International Unit) of a standard hyaluronidase preparation. 1 I-unit = 1 USP-unit = 1 NF-unit.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.