

---

<b>Hygromycin B Lösung</b> Hygromycin B solution	<b>Kat.Nr./Cat.No. :</b>	<b>25966</b>
---	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>20</sub> H <sub>37</sub> N <sub>3</sub> O <sub>13</sub>
<b>Molekulargewicht</b> molecular weight		527.53
<b>Aussehen</b> appearance		klare cognacfarbige Lösung in H <sub>2</sub> O clear cognac-coloured solution in H <sub>2</sub> O
<b>Konzentration (mg/ml)</b> concentration		ca. 400
<b>Mikrobiologische Aktivität (U/ml)</b> microbiological assay		min. 400 000
<b>Gehalt (%)</b> assay	HPLC	min. 80.0
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.