
Inosin reinst Inosine research grade, (6-Hydroxypurine riboside)	Kat.Nr./Cat.No. :	26250
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₂ N ₄ O ₅
Molekulargewicht molecular weight		268.23
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	UV/ 249nm, pH 7.0	min. 98.0
Trocknungsverlust (%) loss on drying		max. 1.0
Schwermetalle (ppm) heavy metals	als / as Pb	max. 10
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.