

Insulin recombinant human	Cat.No. : 26360
-------------------------------------	------------------------

Parameter	Method	Specification
Molecular weight		ca. 5800
Assay (%)	calculated on dry basis	95.0 – 105.0
Assay (U/mg)	USP calculated on dry basis	min. 27.5
Loss on drying (%)		max. 10.0
Bacterial endotoxins (IU/mg)		max. 10
Zinc (%)	calculated on dry basis	0.3 – 0.6
Storage (°C)		-15 to -25

The physical/ chemical data corresponds to the requirements of Ph. Eur. and USP.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.