

Jodnitrotetrazoliumchlorid reinst iodonitrotetrazolium chloride research grade	Kat.Nr./Cat.No. :	26840
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{19}H_{13}N_5O_2Cl$
Molekulargewicht molecular weight		505.72
Aussehen appearance		hellgelbes Pulver light yellow powder
Gehalt (%) assay	Titration	min. 99.0
Löslichkeit solubility	0.5% in Methanol/H ₂ O (1:1)	klare gelbliche Lösung clear yellowish solution
A 1cm/ 249nm	1% in H ₂ O	min. 600
Reduktionsspektrum UV spec. after reducing	Dojindo	entspricht corresponds
Sulfatasche (%) sulfated ash		max. 0.5
Reinheit purity	TLC	einheitlich one spot
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.