

<b>Kinetin p.A.</b> analytical grade	<b>Kat.Nr./Cat.No. :</b>	<b>27190</b>
---	--------------------------	--------------

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>10</sub> H <sub>9</sub> N <sub>5</sub> O
<b>Molekulargewicht</b> molecular weight		215.21
<b>Aussehen</b> appearance		hellelfenbeinfarbenes Pulver light ivory-colored powder
<b>Gehalt (%)</b> assay	HPLC	min. 99.5
<b>A 1cm/ 270nm</b>	0.001% in C <sub>2</sub> H <sub>5</sub> OH	min. 0.75
<b>Schmelzpunkt (°C)</b> melting point		265 - 275
<b>Identität</b> identity	IR	entspricht corresponds
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.