

Kongorot reinst Congo Red research grade, C.I. 22120	Kat.Nr./Cat.No. : 27215
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{32}H_{22}N_6O_6S_2 \cdot Na_2$
Molekulargewicht molecular weight		696.7
Aussehen appearance		braun-rotes Pulver reddish-brown powder
λ_{max} (nm)	0.001% in H ₂ O	498 ± 5
A 1cm/ λ_{max}	0.001% in H ₂ O	min. 0.2
Wasser (%) water	KF	max. 10.0
Reinheit purity	TLC	eine Hauptkomponente, Spuren möglich one major spot, traces possible
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.