

**N-Lauroylsarkosin-Na-salz 30% Lösung**  
N-Lauroylsarcosine Na-salt 30% solution

**Kat.Nr./Cat.No. : 27570**

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>15</sub> H <sub>28</sub> NO <sub>3</sub> Na
<b>Molekulargewicht</b> molecular weight		293.4
<b>Aussehen</b> appearance		gelbliche Flüssigkeit yellowish liquid
<b>Aktive Substanz (%)</b> active substance		28.5 – 32.0
<b>pH</b>	3% in H <sub>2</sub> O	7.5 – 8.5
<b>Fettsäuren (%)</b> fatty acids		max. 2.0
<b>Chlorid (%)</b> chloride	als/ as NaCl	max. 0.2
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.